



Endoscopic Mouthpiece

PATENT NO.: 201120151500.7

PATENT NO.: 201130113710.2

PCT PAT. PEND

Ningbo Greatcare Trading Co., Ltd.

Greatcare's mission is to provide quality medical products with superior value to healthcare providers and end users, improving patient care and the quality of peoples' lives.

GreatCare Medical GmbH was just established in Germany Nov. 2008, the main purpose is to provide after-sale services for all European markets customers, and buy high quality European medical equipments to sell in China markets.

We act with uncompromising integrity. All our dealings are honest. We are committed to treating all people, inside and outside the company, with dignity and respect. We give everyone an equal opportunity, regardless of their personal characteristics, preferences or beliefs.

GreatCare has the strong sales and service team. We have earned the reputation as one of the top sales and customer service companies in China.

GreatCare's new plant has 24,000 square meters. New sales office is located in downtown of Ningbo, China, it has 1,000 square meters, and with 2,000 square meters of warehouse.

As a manufacturer and supplier of several thousands of medical products, GreatCare is dedicated to providing the high quality products. In this pursuit, GreatCare has attained ISO 9001, ISO 13485, EN 46000 compliant, the most stringent international standards for producing quality products. As a result, GreatCare's customers benefit from lower costs and reduced product defects.

With seventeen years of a far-reaching development way and a strong entrepreneurial spirit, Greatcare always focusing on the lofty ideas, developing and making progress of health care, would like to wish the world to create the most trusted future for health care.



Greatcare

**China medical
devices manufacturer,
exporter & importer**

Ningbo Greatcare Trading Co., Ltd.

Unit 93, building 12, No.818 Qiming Rd.

Yinzhou, Ningbo, China.

Tel: +86-574-8778 6632 Fax: +86-574-8778 6644

E-mail: info@greatcare.cn

Creating a Better Care

www.greatcare.cn / www.greaturo.cn

GO

Disposable Endoscopic Mouthpiece

Disposable Endoscopic Mouthpiece with Greatcare patent, is new trend for gastro-fiberoptic endoscope. It consists of two parts: mouthpiece and bandage, and just on this humanistic construction features decrease more pain for patients.



SEE THE FOLLOWING FEATURES IN DETAIL:

The construction of Mouthpiece is hygienic and convenient and prevents transmission of infection disease as well as be adaptable for sleep-induced endoscopies.

The unique feature of the extending tongue depressor not only be comfort for the patients but also propitious to examine.

The groves in the middle prevent the endopiece from moving which allows the teeth to hold it tightly.

The snug fit of the endopiece make patients convenient with its distinguish feature with groove on surface and its polymer material-smooth for the tooth holding.

The novel design of ARC shape of the fixed mount of Endoscopic is suited to human face not to fall off.

The fixed mount with two through holes and buckles cooperate with the bandage flawlessly, that can adjust length according each person. What's more, the bandage with adjust holes is made of silicone material which have patient more comfortable.

The whole Endopiece is made of fully biodegradable plastic; perfectly safe in the environment.

Sterilization method: EO

Package: poly bag



PURPLE MOUTHPIECE



BANDAGE

Disposable Endoscopic Mouthpiece

HOW TO USE:

Let the two ends of bandage through the holes of fixed mount, and fix the bandage with buckles, adjust the bandage length according to each person, then have bandage on patient' head. Put the groves of endopiece into mouth, extending tongue depressor presses the tongue. Patient's teeth are on grooves of endopiece. The endoscopic and the gastroscope through groves to nonnasality, esophagus and tomach to examine.



BEIGE MOUTHPIECE



BEIGE MOUTHPIECE

REF. NO.:	GCR100907
POLY BAG	1PC
MID-BOX	50PCS
CARTON	400PCS

50PCS MOUTHPIECE WITH ONE ELASTIC BANDAGE